

## Accuflo Eclipse (Elastomeric pump)

**Take Medication out of refrigerator and allow to come to room temperature. Check medication label for patient's name, medication and dose, expiration date. **Accuflo devices may vary in shape and size based medication and supplier.****

1. Wash hands/ dry with paper towels
2. Clean your work area as instructed
3. Gather supplies:
  - Medicine in Accuflo (at room temperature)
  - 2 saline flushes
  - 1-heparin flush
  - Alcohol wipes



4. Wash your hands \*\*Prepare Heparin & Saline flushes by removing air from syringes , replacing cap on syringes when done.

5. Unwind the tubing from the Accuflo and remove the cap on the Accuflo tubing but **DO NOT LAY CAP DOWN**

6. Unclamp clam on tubing, allow the fluid to run to the end of the tubing.

7.. Put the cap back on the end of IV tubing and clamp tubing.

9. Clean the injection cap on patient's IV catheter with alcohol x 15 seconds and flush IV catheter with Saline.

10. Clean injection cap on patient's IV catheter with alcohol x 15 seconds and attach IV tubing to the injection cap on patient's IV catheter

11. Unclamp the tubing, be sure **ALL** clamps are open.

12. When the balloon has completely collapsed on the cylinder in the center, the medicine is done.

13. Clamp the IV tubing and wash your hands.

14. Disconnect IV tubing from the injection cap on patient's IV catheter.

15. Clean the injection cap on IV catheter with alcohol x 15 seconds and flush with Saline.

16. Clean the injection cap on patient's IV catheter with alcohol x 15 seconds and flush with Heparin flush

17. Clamp tubing on patient's IV catheter, as instructed.

18. Discard empty Accuflo and syringes as instructed.

